

PRODUCT CATALOG



2024

Company Profile

MAKRITE Industries has been in the business of personal safety and protection for decades — since developing our first exceptional face masks in 1984, well before the global pandemic forced everyone, everywhere to pay attention to respiratory health.

Today our products are in high demand in medical, industrial, automotive and consumer environments worldwide. We have a reliable track record of keeping people safe and productive while keeping costs low.

That's why the National Institute for Occupational Health and Safety (NIOSH) stands behind more than 75 of our disposable respirator products, produced using specially designed, high-tech equipment. This guarantees a superior level of quality and consistency with each face mask we produce. In fact, Makrite products have earned similar approvals—including EN149/CE—from European governments and agencies, along with regulatory bodies in countries such as UK UKCA, Australia AS/NZS, Japan MHLW, Korea KOSHA, China LA and Halal.

We're based in Taiwan, with significant satellite offices in the US, China, Korea and Singapore and Thailand. As the largest private brand manufacturer of disposable respirators for the North American market, Makrite supplies major US retailers and a wide range of other well-known companies that value our reliable customer service, dependable quality and competitive pricing.

After four decades of R&D, we welcome opportunities to further develop our products to meet the changing needs of the industries and partners we serve. We look forward to doing the same for you.

Content

Company Profile	01
NIOSH-Approved Respirators	03
CE & UKCA-Approved Respirators	17
AS/NZS-APPROVED RESPIRATORS	31
JAPANESE MHLW-APPROVED RESPIRATORS	37
KOREAN KOSHA-APPROVED RESPIRATORS	41
HALAL-Approved Respirators	45
LA-Approved Respirators	49
Surgical Masks	53
Makrite Medical Series	59
Makrite & huSHE Series	67
2023 Brand New Products	71

Makrite meets all major global approval standards for respiratory protection.



MAKRITE 9500-N95
TC-84A-5411
NIOSH N95



NIOSH-Approved Respirators


The National Institute for Occupational Safety and Health (NIOSH) is a United States federal agency responsible for conducting research, providing guidance, and making recommendations for the prevention of work-related injuries, illnesses, and fatalities.

Respirators certified by NIOSH means those respirators have met the standards set by NIOSH in terms of filtration efficiency, breathing resistance and air valve leakage.

US NIOSH protective mask grading standards


Grade	Filtering Efficiency			Remark
	95%	99%	99.97%	
N-Type	N95	N99	N100	Not resistant to oil
P-Type	P95	P99	P100	oil Proof
R-Type	R95	R99	R100	Resistant to oil

CERTIFICATE INTRODUCTION:



Approval Holder

Makrite Industries, Inc.
11F-5, No. 79, Sec 1, Xin Tai Wu Road, Xizhi Dist., 221 New Taipei City, Taiwan.
Tel: +886-2-2698-2419



This respirator is approved only in the following configuration:

TC-Approval Number	Filter Designation	Respirator	Cautions and Limitations ²
84A-6806	N95	Envi-910F X	ABCJMNOP

1. Protection

N95-Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply

2. Cautions and Limitations

A - Not for use in atmospheres containing less than 19.5% oxygen.
B - Not for use in atmospheres immediately dangerous to life or health.
C - Do not exceed maximum use concentrations established by regulatory standards.
J - Failure to properly use and maintain this product could result in injury or death.
M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations.
N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
O - Refer to User's Instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
P - NIOSH does not evaluate respirators for use as surgical masks.

NIOSH-APPROVED RESPIRATORS

9600 COMFORT SERIES

The 9600 Comfort Series, our newest line of NIOSH-approved respirators ensure that you are getting a respirator that meets the highest standards of performance and reliability. The foldable design of the Makrite 9600 brings a new level of versatility to your workplace. Easily stow it away when not in use, and effortlessly unfold it when you need immediate protection. Its compact size makes it ideal for professionals on the go, ensuring you're always prepared for any situation.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
9600-N95	N95 Cone Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-7867
9600-N95S	N95 Cone Respirator	S	20 masks/box; 12 boxes/case	TC-84A-7871
9600V-N95	N95 Cone Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-7869
9600V-N95S	N95 Cone Respirator with exhalation valve	S	10 masks/box; 12 boxes/case	TC-84A-8107
9600-N95OV	N95 Cone Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	TC-84A-7882
9600-N95VOV	N95 Cone Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	TC-84A-7883
9601-N95	N95 Cone Respirator with adjustable buckle	M/L	20 masks/box; 12 boxes/case	TC-84A-9225
9601-N95S	N95 Cone Respirator with adjustable buckle	S	20 masks/box; 12 boxes/case	TC-84A-9287
9606	N95 Cone Respirator with black outer layer	M/L	20 masks/box; 12 boxes/case	TC-84A-9463
9606V	N95 Cone Respirator w/vlave and black outer layer	M/L	10 masks/box; 12 boxes/case	TC-84A-9498

NIOSH-APPROVED RESPIRATORS

9500 COMFORT SERIES

Experience unmatched protection with the Makrite 9500, the leading series of NIOSH-approved N95 disposable respirators from Makrite. Trusted by professionals worldwide, the Makrite 9500 offers superior respiratory protection in a comfortable cone style design. With its high filtration efficiency, it effectively blocks out airborne particles, ensuring clean and safe breathing. The adjustable nose piece and secure fit provide added comfort and confidence throughout the day.




ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
* 9500-N95	N95 Cone Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-5411
* 9500-N95S	N95 Cone Respirator	S	20 masks/box; 12 boxes/case	TC-84A-5463
9500V-N95	N95 Cone Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-5460
9500-N95OV	N95 Cone Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	TC-84A-5525
9500-N95VOV	N95 Cone Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	TC-84A-6954
8810S	N95 Cone Respirator (looping style headstrap)	S	20 masks/box; 12 boxes/case	TC-84A-6953
9500I	N95 Cone Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-9420
9500SI	N95 Cone Respirator	S	20 masks/box; 12 boxes/case	TC-84A-9451

* Cleared by the Food and Drug Administration for use as a Surgical N95 Respirator, FDA 510(k): K020474

NIOSH-APPROVED RESPIRATORS

MK COMFORT SERIES

A premium choice for those seeking an alternative to our 9600 and 9500 models. The MK Comfort series features a distinctive flanged edge encircling the mask, ensuring a snug and secure fit for some users. In response to wearer preferences, we have reimagined the internal nosepiece, crafting it with a thinner, elongated strip of PVC foam. Setting it apart from our 9600 and 9500 Comfort Series, the MK Comfort Series is equipped with a square exhalation valve, providing an added element of differentiation. Elevate your comfort and confidence with the MK Comfort Series, setting new standards of excellence in respiratory protection.



MK910-N95



MK910-N95S



MKN95-910V

NIOSH-APPROVED RESPIRATORS

ENVI COMFORT SERIES

A highly sought-after choice in the food service industry renowned for its eco-friendly composition, free from any metallic components. The Envi Comfort Series features a soft, pre-shaped nose area that guarantees an exceptional fit across diverse face sizes and contours, effortlessly conforming to the wearer's unique nasal bridge for optimal comfort. For added convenience, the looping-style headstrap enable easy donning and doffing of the mask, while also allowing it to be conveniently draped around the neck when not in use. Complemented by a soft inner layer and optimized breathing resistance, the Envi Series offers unparalleled comfort that surpasses expectations.



Envi910



Envi910V



Envi-910F



Envi-910FV



N95-710VOV

ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
MK910-N95	N95 Cone Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-5408
MK910-N95S	N95 Cone Respirator	S	20 masks/box; 12 boxes/case	TC-84A-5464
MKN95-910V	N95 Cone Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-5409
MK910-N95OV	N95 Cone Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	TC-84A-7829
MK910-N95VOV	N95 Cone Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	TC-84A-7830

ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
Envi910	N95 Cone Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-4656
Envi910V	N95 Cone Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-4657
Envi-910F	N95 Cone Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-6806
Envi-910FV	N95 Cone Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-6811
N95-710VOV	N95 Cone Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	TC-84A-9219

NIOSH-APPROVED RESPIRATORS

VERTICAL FOLD SERIES

In North and South America,, there is a noticeable trend among wearers, with approximately 25% favoring the versatility and adaptability of vertical- and flat-folded respirators—a figure that continues to rise annually. Many individuals prefer these styles due to their enhanced flexibility around the edges, allowing for a seamless fit that conforms to facial features with greater ease compared to cone-style masks. When not in use, this mask conveniently folds down the middle, facilitating effortless storage. Furthermore, it fits well with goggles or helmets, making it an ideal choice for a wide range of work environments.



SEKURA-N95 SEKURA-511 MAK 368+ MAK 699+

ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
SEKURA-321	N95 Vertical Folded Respirator with adjustable buckle	M/L	20 masks/box; 12 boxes/case	TC-84A-6660
SEKURA-511	N95 Vertical Folded Respirator w/valve and adjustable buckle	M/L	10 masks/box; 12 boxes/case	TC-84A-6666
SEKURA-N95	N95 Vertical Folded Respirator	M/L	40 masks/box; 12 boxes/case	TC-84A-7638
SEKURA-N95V	N95 Vertical Folded Respirator with exhalation valve	M/L	15 masks/box; 12 boxes/case	TC-84A-7643
MAK368 ⁺	N95 Vertical Folded Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-7774
MAK520 ⁺	N95 Vertical Folded Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-7773
MAK699 ⁺	N95 Vertical Folded Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-7841

NIOSH-APPROVED RESPIRATORS

VERTICAL FOLD SERIES & FLAT FOLD SERIES

In North and South America, there is a discernible trend wherein approximately 25% of consumers favor the vertical and flat folded style respirator, a figure that continues to experience steady growth. This preference is driven by the respirator's exceptional flexibility along the edges, enabling seamless adaptation to the wearer's facial features. In fact, for some individuals, the vertical and flat folded style offers even greater adaptability compared to traditional cone-style respirators.



910-N95FMX N95V-910FM 6001 801-N95 2201

ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
* 910-N95FMX	N95 Vertical Folded Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-4541
N95V-910FM	N95 Vertical Folded Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-4542
6001	N95 Vertical Folded Respirator	M/L	30 masks/box; 12 boxes/case	TC-84A-9372
6001V	N95 Vertical Folded Respirator with exhalation valve	M/L	15 masks/box; 12 boxes/case	TC-84A-9435
801-N95	N95 Vertical Folded Respirator w/adjustable buckle & M-shaped foam	M/L	10 masks/box; 12 boxes/case	TC-84A-9223
2201	N95 Vertical Folded Respirator	M/L	30 masks/box; 12 boxes/case	TC-84A-9230
2201V	N95 Vertical Folded Respirator with exhalation valve	M/L	20 masks/box; 12 boxes/case	TC-84A-9231

* Cleared by the Food and Drug Administration for use as a Surgical N95 Respirator, FDA 510(k): K020474

NIOSH-APPROVED RESPIRATORS

500 COMFORT SERIES

Specifically designed to protect against metal fumes, including ozone. These exceptional respirators feature a very smooth inner surface, ensuring excellent comfort during extended wear. Equipped with fully adjustable headstraps and a full inner face seal, they guarantee a secure and personalized fit for every wearer. The flame retardant outer layer adds an extra layer of protection, while the preformed nose area, enhanced with an adjustable aluminum nose clip, ensures a seamless fit and optimal wearer comfort.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
510-N95W	N95 Welding Respirator with adjustable four-point headstraps	M/L	10 masks/box; 12 boxes/case	TC-84A-4654
N95V-510W	N95 Welding Respirator w/valve and adjustable four-point headstraps	M/L	10 masks/box; 12 boxes/case	TC-84A-4652
N95VOV-510W	N95 Welding Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	TC-84A-6955
MKN95V-530W	N95 Welding Respirator w/valve and headstraps	M/L	10 boxes/case;12 boxes/case	TC-84A-8074
MKN95V-530Wa	N95 Welding Respirator w/valve and adjustable four-point headstraps	M/L	10 boxes/case;12 boxes/case	TC-84A-8076

NIOSH-APPROVED RESPIRATORS

N99 RESPIRATORS

The N99 provides greater than 99% filter efficiency protection against non-oil based particulates. Our NIOSH approved N99 respirators contain breathable, skin-friendly and soft properties, which can be used in various settings and works against air pollutants, pollen and dust particles.




ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
N99-810V	N99 Respirator with exhalation valve(stapled headstraps)	M/L	10 masks/box; 12 boxes/case	TC-84A-4081
N99-820V	N99 Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-7639

NIOSH-APPROVED RESPIRATORS

P95 RESPIRATORS

P95 respirators provide greater than 95% filter efficiency against oil and non-oil-based particulate aerosols. The P95 is strongly resistant to oils, according to NIOSH standards. These respirators are also available with a layer of activated carbon for nuisance odor relief and an exhalation valve for easy breathing.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
P95-710VOV	P95 Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	TC-84A-5041
MK710-P95VOV	P95 Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	TC-84A-7528
MK710-P95V	P95 Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	TC-84A-7696

NIOSH-APPROVED R95 RESPIRATORS

R95 RESPIRATORS

R95 respirators provide greater than 95% filter efficiency against oil and non-oil-based particulate aerosols. The R95 is somewhat resistant to oils, according to NIOSH standards. Different versions are also available with a layer of activated carbon for nuisance odor relief and an exhalation valve for easy breathing.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
R95-610V	R95 Respirator with valve and full PVC foam	M/L	10 masks/box; 12 boxes/case	TC-84A-5154
R95-620V	R95 Respirator with valve and half PVC foam	M/L	10 masks/box; 12 boxes/case	TC-84A-5437
MK610-R95OV	R95 Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	TC-84A-7619

NIOSH-APPROVED N95 RESPIRATORS

HYBIRD STYLE

Our innovative FIT-N95 respirator, available with or without an exhalation valve, combines the advantages of cone and flat folded styles in one. Engineered with flexible edges, this respirator conforms to the facial contours of most wearers, ensuring a secure yet comfortable seal. The unique fold-out front design provides distance from the mouth, similar to traditional cone-style masks, enhancing breathability and reducing moisture buildup.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
FIT-N95	N95 Respirator	M/L	30 masks/case;12 boxes/case	TC-84A-8420
FIT-N95V	N95 Respirator with exhalation valve	M/L	15 masks/case;12 boxes/case	TC-84A-8421

NIOSH-APPROVED N95 RESPIRATORS

9800 SERIES

Our commitment to environmental protection is evident in every aspect of this innovative design. A reusable full face seal series with food-grade silicone material was developed for this purpose. The soft silicone gasket integrated into each respirator enhances both comfort and fit, ensuring a secure and customized seal for every wearer.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
9800-N95	N95 Respirator	M/L	20 masks/case;12 boxes/case	TC-84A-9049
9800V-N95	N95 Respirator with exhalation valve	M/L	10 masks/case;12 boxes/case	TC-84A-9055

CE & UKCA-APPROVED RESPIRATORS

The CE certification is an essential requirement for products intended to be sold in the European Economic Area (EEA). It serves as a mark of compliance with rigorous European health, safety, and environmental protection standards. When you see the CE logo on a commercial product, it signifies that the manufacturer or importer affirms the product's conformity to these standards. In the case of masks, which fall under the category of personal protective equipment (PPE), compliance with the Personal Protective Equipment Regulation (EU) 2016/425 is crucial. Rest assured, Makrite respirators meet the stringent CE standards, ensuring their ability to effectively filter both non-oil and oil-based particles for optimal protection.

FFP1- filters at least 80% of airborne particles
FFP2- filters at least 94% of airborne particles
FFP3- filters at least 99% of airborne particles


UKCA is the abbreviation of UK Conformity Assessed. The UKCA mark is a safety mark for British products and is applicable to products sold in England, Scotland and Wales.



CE & UKCA-APPROVED RESPIRATORS

STANDARD COMFORT SERIES

The Standard Comfort Series is officially approved for use in European Union countries and others that follow the CE respirator protocols. These respirators meet all levels of EU safety standards, including FFP1, FFP2, and FFP3 classifications, ensuring optimal protection for the wearer. They are compatible with various head- and hearing-protection products, offering versatility in different work environments. Our range includes models without exhalation valves, suitable for work environments that pose no breathing challenge, as well as options with advanced exhalation valves for enhanced comfort in hot and humid work settings. Additional features such as adjustable headstraps, various comfort enhancements, and an optional carbon layer for odor relief are available, allowing for a tailored and personalized experience.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
FP1SLw	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP1SLpw	FFP1 NR Respirator (PP outer layer)	M/L	20 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP1SLAVw	FFP1 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP1SLAVpw	FFP1 NR Respirator with exhalation valve (PP outer layer)	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP1SLOVw	FFP1 NR Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP1SLOVpw	FFP1 NR Respirator for nuisance organic vapors (PP outer layer)	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP1VSLOVw	FFP1 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP1VSLOVpw	FFP1 NR Respirator for nuisance organic vapors w/valve (PP outer layer)	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP2SLw	FFP2 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP2SLpw	FFP2 NR Respirator (PP outer layer)	M/L	20 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP2SLAVw	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP2SLAVpw	FFP2 NR Respirator with exhalation valve (PP outer layer)	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP2SLOVw	FFP2 NR Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP2SLOVpw	FFP2 NR Respirator for nuisance organic vapors (PP outer layer)	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP2VSLOVw	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP2VSLOVpw	FFP2 NR Respirator for nuisance organic vapors w/valve (PP outer layer)	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP3SL	FFP3 NR Respirator with stapled headstrap	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP3SLAV	FFP3 NR Respirator w/valve and adjustable straps	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP3SLw	FFP3 NR Respirator	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
FP3SLAVpw	FFP3 NR Respirator with exhalation valve (PP outer layer)	M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308
		M/L	10 masks/box; 12 boxes/case	CE 711366 & UKCA759308

CE & UKCA-APPROVED RESPIRATORS

MK COMFORT SERIES

Our MK Comfort Series, an exceptional alternative to our standard/original CE cone respirator. Building upon the same outstanding features and advantages of the original model, this series incorporates a flanged edge around the entire mask. This innovative design ensures a secure and tight seal with the face, catering to the preferences of discerning wearers seeking an optimal fit. Experience the comfort and performance excellence of our MK Comfort Series respirators.




ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
MKP1SLw	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 712994 & UKCA 759311
MKP1SLAVw	FFP1 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712994 & UKCA 759311
MKP2SLa	FFP2 NR Respirator with adjustable straps	M/L	20 masks/box; 12 boxes/case	CE 712994 & UKCA 759311
MKP2SLAVa	FFP2 NR Respirator w/valve and adjustable straps	M/L	10 masks/box; 12 boxes/case	CE 712994 & UKCA 759311
MKP2SLw	FFP2 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 712994 & UKCA 759311
MK2SLAVw	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712994 & UKCA 759311
MKP2VSLOVw	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 712994 & UKCA 759311
MKP3SLAVa	FFP3 NR Respirator w/valve and adjustable straps	M/L	10 masks/box; 12 boxes/case	CE 712994 & UKCA 759311
MKD110	FFP1 NR D Respirator	M/L	20 masks/box; 12 boxes/case	CE 746651 & UKCA 759313
MKD110V	FFP1 NR D Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 746651 & UKCA 759313
MKD220	FFP2 NR D Respirator	M/L	20 masks/box; 12 boxes/case	CE 746651 & UKCA 759313
MKD220V	FFP2 NR D Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 746651 & UKCA 759313
MKD220Va	FFP2 NR D Respirator w/valve and adjustable straps	M/L	10 masks/box; 12 boxes/case	CE 746651 & UKCA 759313
MKD330Va	FFP3 NR D Respirator w/valve and adjustable straps	M/L	10 masks/box; 12 boxes/case	CE 746651 & UKCA 759313
MKD330Vra	FFP3 NR D Respirator w/valve and adjustable straps	M/L	10 masks/box; 12 boxes/case	CE 746651 & UKCA 759313

CE & UKCA-APPROVED RESPIRATORS

FLAT FOLD SERIES

Particularly popular in Europe, the flat-fold style respirator offers a flexibility in fit that many wearers prefer. Experience the exceptional of comfort and protection as the soft and flexible edges seamlessly conform to the unique facial features of each individual, ensuring an optimal fit that enhances both comfort and performance.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
FP1SLFM	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP1VSLFMc	FFP1 NR Respirator with exhalation valve (concealed nose clip)	M/L	20 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP1VSLFM	FFP1 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP1VSLOVFM	FFP1 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP2SLFM	FFP2 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP2VSLFM	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP2SLOVFM	FFP2 NR Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP2VSLOVFM	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP2VSLOVFMc	FFP2 NR Respirator for nuisance organic vapors w/valve (concealed nose clip))	M/L	10 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP3SLFM	FFP3 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 712989 & UKCA 759309
FP3VSLFM	FFP3 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712989 & UKCA 759309

CE & UKCA-APPROVED RESPIRATORS

CFP SERIES

Our CFP Comfort Series, thoughtfully designed with foldable respirators that effortlessly fit into a worker's pocket, ensuring convenient storage and accessibility. The unique shape of these respirators ensures an optimal fit for a diverse range of individuals in any workplace, providing enhanced comfort and performance. These models exceed the rigorous CE approval standards for breathing resistance and filtration efficiency, offering the utmost reliability and protection




ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
CFP1SL	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 711366&UKCA 759308
CFP2SL	FFP2 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 711366&UKCA 759308
CFP2SLAV	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 711366&UKCA 759308
CFP2SLOV	FFP2 NR Respirator for nuisance organic vapors	M/L	20 masks/box; 12 boxes/case	CE 711366&UKCA 759308
CFP2SLVOV	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 711366&UKCA 759308
CFP3SLAV	FFP3 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 711366&UKCA 759308

CE & UKCA-APPROVED RESPIRATORS

9800 SERIES

Our revolutionary line of reusable respirators, designed with a strong commitment to environmental preservation and minimizing product-related pollution. The 9800 series features a food-grade silicone full face seal, enhanced with a soft edge for superior protection and all-day wearer comfort.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
9800P1	FFP1 NR Respirator	M/L	20 boxes/case;12 boxes/case	CE 712994 & UKCA 7559311
9800P1V	FFP1 NR Respirator with exhalation valve	M/L	10 boxes/case;12 boxes/case	CE 712994 & UKCA 7559311
9800P2	FFP2 NR Respirator	M/L	20 boxes/case;12 boxes/case	CE 712994 & UKCA 7559311
9800P2V	FFP2 NR Respirator with exhalation valve	M/L	10 boxes/case;12 boxes/case	CE 712994 & UKCA 7559311
9800P2OV	FFP2 NR Respirator for nuisance organic vapors	M/L	10 boxes/case;12 boxes/case	CE 712994 & UKCA 7559311
9800P2VOV	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	10 boxes/case;12 boxes/case	CE 712994 & UKCA 7559311

CE & UKCA-APPROVED RESPIRATORS

FIT SERIES

Our FIT Comfort Series (available with or without an exhalation valve) combines the best elements of cone-style and flat-folded respirators. Designed with precision, the flexible edges effortlessly conform to the wearer's facial contours, ensuring a secure yet comfortable fit. Moreover, the fold-out front provides the desired distance from the mouth just as the traditional cone-style respirators.




ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
FIT-P1DSL	FFP1 NR D Respirator	M/L	30 bags/box;12 boxes/case	CE 763157 & UKCA 776079
FIT-P1DSLAV	FFP1 NR D Respirator with exhalation valve	M/L	15 bags/box;12 boxes/case	CE 763157 & UKCA 776079
FIT-P2DSL	FFP2 NR D Respirator	M/L	30 bags/box;12 boxes/case	CE 763157 & UKCA 776079
FIT-P2DSLAV	FFP2 NR D Respirator with exhalation valve	M/L	15 bags/box;12 boxes/case	CE 763157 & UKCA 776079
FIT-P2DSLVOV	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	15 bags/box;12 boxes/case	CE 763157 & UKCA 776079
FIT-P3DSL	FFP3 NR D Respirator	M/L	30 bags/box;12 boxes/case	CE 763157 & UKCA 776079
FIT-P3DSLAV	FFP3 NR D Respirator with exhalation valve	M/L	15 bags/box;12 boxes/case	CE 763157 & UKCA 776079
FIT-P3DSLAVa	FFP3 NR D Respirator with exhalation valve	M/L	15 bags/box;12 boxes/case	CE 763157 & UKCA 776079

CE & UKCA-APPROVED RESPIRATORS

500 COMFORT SERIES

The 500 Comfort Series is designed with a prominent emphasis on environmental friendliness and comfort. Its unique pre-shaped nose area guarantees an exceptional fit for diverse facial sizes and shapes. Free from any metallic components, the series incorporates a soft inner layer, offering users the utmost comfort during wear. Available in both FFP1 and FFP2, it also provides options with or without an exhalation valve and an additional activated carbon layer, meeting the specific work environments and individual requirements.




ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
Envi 501P1	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
Envi 501P1OV	FFP1 NR Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
Envi 501P1V	FFP1 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
Envi 501P1 VOV	FFP1 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
Envi 502P2	FFP2 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
Envi 502P2OV	FFP2 NR Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
Envi 502P2V	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
Envi 502P2VOV	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312

CE & UKCA-APPROVED RESPIRATORS

ENVI COMFORT SERIES

The ENVI Series of respirators is designed without any metal in the construction, ensuring they are environmentally friendly. The soft inner layer provides maximum comfort for the user's face. The adjustable looping-style head straps offer a convenient way for individual adjustment during donning and doffing. The head strap loop design allows the mask to be easily hung around the neck when not in use. With low breathing resistance, these respirators offer a highly comfortable wear experience. The soft pre- shaped nose area ensures an excellent fit for various facial sizes and should be adjusted up and down the wearer's nose bridge area for a proper fit. Due to these features, this respirator has gained significant preference in the food service industry.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
ENVIP1F	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP2VF	FFP2 NR R Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP3VF	FFP3 NR R Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP1DF	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP1DVF	FFP1 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP2DF	FFP2 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP2DVF	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP2DVOVF	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP3DVF	FFP3 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP3DF	FFP3 NR Respirator	M/L	10 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP3Da+	FFP3 NR D Respirator with adjustable strap	L	20 masks/box; 12 boxes/case	CE 712990 & UKCA 759310
ENVIP3DVa+	FFP3 NR D Respirator w/valve and adjustable strap	L	10 masks/box; 12 boxes/case	CE 712990 & UKCA 759310

CE & UKCA-APPROVED RESPIRATORS

300 COMFORT SERIES

Another version of the foldable respirator is the vertical fold. These lightweight respirators offer clear vision with no possibility of obstructed views. They are available in a wide range of versions as shown in the chart below, providing protection for a variety of industrial situations.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
6001-P2	FFP2 NR Respirator	M/L	40 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
301P1	FFP1 NR Respirator	M/L	20 boxes/case;12 boxes/case	CE 718846 & UKCA 759312
301P1V	FFP1 NR Respirator with exhalation valve	M/L	10 boxes/case;12 boxes/case	CE 718846 & UKCA 759312
301P1w	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
301P1OVw	FFP1 NR Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
301P1Vw	FFP1 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
301P1VOVw	FFP1 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
302P2	FFP2 NR Respirator	M/L	20 boxes/case;12 boxes/case	CE 718846 & UKCA 759312
302P2V	FFP2 NR Respirator with exhalation valve	M/L	10 boxes/case;12 boxes/case	CE 718846 & UKCA 759312
302P2w	FFP2 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
302P2OVw	FFP2 NR Respirator for nuisance organic vapors	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
302P2Vw	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
302P2VOVw	FFP2 NR Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
303P3w	FFP3 NR Respirator	M/L	20 masks/box; 12 boxes/case	CE 718846 & UKCA 759312
303P3Vw	FFP3 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	CE 718846 & UKCA 759312

CE & UKCA-APPROVED RESPIRATORS

300D & 880 COMFORT SERIES

Compare to the 300 comfort series, the head-straps of 300D series masks are replaced with anti-skidding rubber material, making it more convenient for the wearers to adjust the straps with the adjustable buckles design. Additionally, the mask features fully adjustable loop head-straps for proper fit. The new 880 Comfort Foldable Mask Series meets the optional requirement of dolomite clogging test. The mask design improves the fit between the mask and the face, providing enhanced safety protection. In addition, the streamlined patterns on the mask provide a modern aesthetic. The foldable design allows for easy storage and portability and can effectively be used in conjunction with safety goggles and helmets.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
301P1Da	FFP1 NR D Respirator (dolomite tested)	M/L	20 masks/box; 12 boxes/case	190974A
302P2DVa	FFP2 NR D Respirator with exhalation valve (dolomite tested)	M/L	10 boxes/case;12 boxes/case	190975A
303P3DVa	FFP3 NR D Respirator with exhalation valve (dolomite tested)	M/L	10 boxes/case;12 boxes/case	190976A
303P3VWPLUS	FFP3 NR Respirator w/valve and adjustable buckle	M/L	10 boxes/case;12 boxes/case	190502A
880-FFP1D	FFP1 NR D Respirator (dolomite tested)	M/L	30 boxes/case;12 boxes/case	CE 718846 & UKCA 759312
880-FFP2D	FFP2 NR D Respirator (dolomite tested)	M/L	30 boxes/case;12 boxes/case	CE 718846 & UKCA 759312
880-FFP2DV	FFP2 NR D Respirator with exhalation valve (dolomite tested)	M/L	15 boxes/case;12 boxes/case	CE 718846 & UKCA 759312
880-FFP3DV	FFP3 NR D Respirator with exhalation valve (dolomite tested)	M/L	15 boxes/case;12 boxes/case	CE 718846 & UKCA 759312



CE & UKCA-APPROVED RESPIRATORS

ET SERIES

The ET Series is the economical version of our Standard Comfort Series, designed to meet highly competitive pricing demands. Despite its straightforward design, its performances and applications are no different to our other CE and UKCA approved respirators.

individually adjustable
nose clip for excellent fit

ultrasonically welded
headstrap for secure seal

low breathing resistance
for comfort



FFP1

FFP1

FFP2

FFP2

FFP2

FFP2



ET-X P1VN

ET-X P1OVN

ET-X P2N

ET-X P2VN

ET-X P2OVN

ET-X P2VOVN

ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
ET-X P1N	FFP1 NR Respirator	M/L	30 masks/box; 12 boxes/case	CE 711366 & UKCA 759308
ET-X P1VN	FFP1 NR Respirator with exhalation valve	M/L	15 boxes/case;12 boxes/case	CE 711366 & UKCA 759308
ET-X P1OVN	FFP1 NR for nuisance organic vapors	M/L	20 masks/box; 12 boxes/case	CE 711366 & UKCA 759308
ET-X P2N	FFP2 NR Respirator	M/L	30 masks/box; 12 boxes/case	CE 711366 & UKCA 759308
ET-X P2VN	FFP2 NR Respirator with exhalation valve	M/L	15 boxes/case;12 boxes/case	CE 711366 & UKCA 759308
ET-X P2OVN	FFP2 NR for nuisance organic vapors	M/L	20 masks/box; 12 boxes/case	CE 711366 & UKCA 759308
ET-X P2VOVN	FFP2 NR for nuisance organic vapors w/valve	M/L	15 boxes/case;12 boxes/case	CE 711366 & UKCA 759308



AS/NZS-APPROVED RESPIRATORS

In the Australia and New Zealand regions, this standard is required to be complied with by workers who are using respiratory protective equipment. The standard specifies performance requirements and testing criteria for respiratory protective devices intended to provide protection against atmospheres containing substances that may be harmful if breathed.

bsi. 
By Royal Charter

Certificate

This is to certify that **Makrite Industries Inc.**

ABN Not Applicable
Head Office Address 12F-5, No. 79, Sec. 1, Xintai 5th Road, Xizhi Dist., New Taipei City 221, Taiwan
Manufacturing Site Makrite Industries Inc. DONGGUAN SURT, ON PLASTIC & FIBER PRODUCTS CO., LTD. Zhufu Ling Industrial Area, Tangsha Town, Dongguan City, Guangdong, CHINA.

Holds Certificate Number **BMP 720690**

Is authorized to use the BSI registered Benchmark Product Certification Trademark, on products identified in the Product Certification Schedule that is attached to this Product Certification Certificate. Such products have been Type Tested and are manufactured under the Product Certification Requirements monitored by BSI to ensure that the manufacturing process has the capability to consistently produce products in compliance with and are Certified to:

AS/NZS 1716:2012 - Respiratory protective devices

For and on behalf of BSI 
Hao Barnes, Managing Director, BSI Group ANZ Pty Ltd

Originally registered 24/02/2009 Latest issue 02/03/2023 Expiry date 01/03/2025


  **...making excellence a habit**

The use of this certificate and of the Benchmark Product Certification Trademark is subject to the terms and conditions of the BSI Benchmarking and Certification Scheme. BSI Group ANZ Pty Ltd, Suite 1, 14 Robinson Road, Singapore 048533. The validity of this certificate can be checked at www.bsi.com.au/certification and the BSI website www.bsi.com.au. This certificate remains the property of BSI and is not to be reproduced without the written permission of BSI Group ANZ Pty Ltd.

AS/NZS-APPROVED RESPIRATORS

MK STANDARD COMFORT SERIES

The MK Standard Comfort Series features a soft inner layer, an adjustable nose clip, and a flanged edge, ensuring a secure and reliable facial seal. With durability in mind, the ultrasonically welded headstraps offer a strong seal that can withstand heavy use.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
MKP1SLw	P1 Respirator	M/L	20 masks/box; 12 boxes/case	BMP 720690
MKP1SLAVw	P1 Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	BMP 720690
MKP2SLw	P2 Respirator	M/L	20 masks/box; 12 boxes/case	BMP 720690
MKP2SLAVw	P2 Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	BMP 720690
MKP2VSLOVw	P2 Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	BMP 720690

AS/NZS-APPROVED RESPIRATORS

FIT & 300 COMFORT SERIES

Our FIT Comfort Series (available with or without an exhalation valve) combines features of both the cone and flat-folded styles of respirators. The flexible edges conform to wearers’ facial features with a tight but comfortable seal, and the fold-out front allows for the same distance from the mouth as in traditional cone-style masks. Soft to the touch, the 300 Comfort Series is exceptionally easy to wear. Ultrasonically welded headstraps provide a strong seal and can withstand heavy use.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
FIT-P2SL	P2 Respirator	M/L	20 masks/box; 12 boxes/case	BMP 720690
FIT-P2SLAV	P2 Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	BMP 720690
301P1w	P1 Respirator	M/L	30 masks/box; 12 boxes/case	BMP 720690
302P2Vw	P2 Respirator with exhalation valve	M/L	15 masks/box; 12 boxes/case	BMP 720690

AS/NZS-APPROVED RESPIRATORS

Envi COMFORT SERIES

Especially popular in the food service industry, the Envi Series is made without the use of any metals, making it environmentally friendly. The soft, pre-shaped nose area ensures excellent fit for various face sizes and shapes, adjusting to the wearer’s nose bridge for a proper fit. In addition, the mask’s soft inner layer and low breathing resistance ensure maximum comfort.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
Envi 502P2	P2 Respirator	M/L	20 masks/box; 12 boxes/case	BMP 720690
Envi 502P2V	P2 Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	BMP 720690
Envi 502P2VOV	P2 Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	BMP 720690

AS/NZS-APPROVED RESPIRATORS

ET COMFORT SERIES

As the value-based line of P1 and P2 models, the ET Comfort Series is designed to meet very competitive pricing demands. Their function and uses are the same as our other AS/NZS-approved respirators.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
ET-P1	P1 Respirator	M/L	20 masks/box; 12 boxes/case	BMP 720690
ET-P1VOV	P1 Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	BMP 720690
ET-P2	P2 Respirator	M/L	20 masks/box; 12 boxes/case	BMP 720690
ET-P2V	P2 Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	BMP 720690
ET-P2OV	P2 Respirator for nuisance organic vapors	M/L	20 masks/box; 12 boxes/case	BMP 720690
ET-P2VOV	P2 Respirator for nuisance organic vapors w/valve	M/L	10 masks/box; 12 boxes/case	BMP 720690



JAPANESE MHLW-APPROVED RESPIRATORS


The Japanese MHLW (Ministry of Health, Labour and Welfare) is a comprehensive organization that oversees the general health and welfare of citizens, which includes occupational safety and health. As such, it governs the standards for occupational health and safety equipment, including respirators. It is the official authority responsible for issuing approvals for Japanese standard products.

労働衛生保護具型式検定合格証		産業安全技術協会 JSA Japanese Society of Safety Engineering	
申請者	11F-5, No.79, Sec.1, XinTai Wu Rd., XIZHI Dist. New Taipei city, Taiwan, R.O.C. Makino Industries Inc.		
製造者	Building 2, Zhuo Feng Ling 168, Industries Area, TangXia, Dongguan City, Guangdong, China Saw-ven Factory, Makino Industries Inc.		
品名	防じんマスク		
型式の名称	MK910-N95DS2		
型式検定合格番号	第 TM531 号	区	分 DS2
構造	顔当て式 (平面形)		
使用条件			
有効期間		型式検定者の所属及び氏名	
2017年 7月12日から 2022年 7月11日まで		試験認証部 熊島 直之	
2022年 7月12日から 2027年 7月11日まで		試験認証部 熊島 直之	
年 月 日から	年 月 日まで		
年 月 日から	年 月 日まで		
機械等検定規則による型式検定に合格したことを証明する。 2012年 7月12日			
型式検定実施者 公益社団法人 産業安全技術協会会長			
TM531_8.1/1			

MHLW-APPROVED RESPIRATORS

MHLW SERIES

The MK 910-DS2 series respirators are certified by the Japanese MHLW standard. These respirators are designed as an alternative to our 9500 comfort series. The primary difference between the two is that the MK design features a "flanged" edge all around the mask, as some users prefer this type of seal to the face. Additionally, the MK respirator includes a PVC nose foam on the inner surface, providing a better seal. The filter media is designed for low breathing resistance, ensuring maximum comfort for users.



ITEM	SIZE	PACKAGING	APPROVAL
MK910-N95DS2	M/L	20 masks/box; 12 boxes/case	TM531
MK910-N95DS2V	M/L	10 masks/box; 12 boxes/case	TM532

MHLW-APPROVED RESPIRATORS

MHLW SERIES

The 9600 DS2 Series, our newest collection of MHLW-approved respirators ensure that you are getting a respirator that meets the highest standards of performance and reliability. The foldable design of the 9600 DS2 brings a new level of versatility to your workplace. Easily stow it away when not in use, and effortlessly unfold it when you need immediate protection. Its compact size makes it ideal for professionals on the go, ensuring you're always prepared for any situation.




ITEM	SIZE	PACKAGING	APPROVAL
9600-DS2	M/L	20 masks/box; 12 boxes/case	TM814
9600-DS2V	M/L	10 masks/box; 12 boxes/case	TM815



KOREAN KOSHA-APPROVED RESPIRATORS

The Korea Occupational and Health Agency was founded with an aim to build a safer and healthier working environment for all. KOSHA strives to prevent occupational accidents and diseases, promote workers safety and health, and encourage enterprises to launch accident prevention activities.





KOSHA-APPROVED RESPIRATORS

KOSHA SERIES

We have a range of products that comply with the Korean national standards. These products range from three different efficiency levels. As we aim to provide more diversity to the market, we have chosen various shapes and designs to cover the needs of a wider range of users.

low breathing resistance for comfort

adjustable concealed plastic nose clip for proper fit

one-way valve for easy exhalation



SLAV-301



301V-103



501-102



SLw-101



SLAVw-201



VSLOVw-202

ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
301V-103	FFP1 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	11-AVACT-0123
501-102	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	11-AV4CT-0124
SLw-101	FFP1 NR Respirator	M/L	20 masks/box; 12 boxes/case	11-AV4CT-0125
SLAVw-201	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	11-AV4CT-0120
VSLOVw-202	FFP2 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	11-AV4CT-0121
SLAVw-301	FFP3 NR Respirator with exhalation valve	M/L	10 masks/box; 12 boxes/case	11-AV4CT-0119

HALAL-Approved Respirators

"Halal" is an Arabic word that translates to 'permissible' in English. A Halal certificate is a qualification issued by a qualified independent third-party supervisory organization for food, medicine, chemicals, additives, etc., ensuring that they conform to Halal regulations as required by Muslims. A product with a Halal certificate indicates that it is suitable and acceptable for Muslims.



HALAL-APPROVED RESPIRATORS

Find more details on the referring page



9600-N95 (P5)



9600-N95S (P5)



9600V-N95 (P5)



9500-N95 (P6)



9500-N95S (P6)



9500V-N95S (P6)



SEKURA-N95 (P9)



SEKURA-N95V (P9)



SEKURA-321 (P9)



SEKURA-511 (P9)



MAK368+ (P9)



E101W (P65)



E103B (P65)



P603BE (P66)



LP603BE (P66)





LA-Approved Respirators

In 2015, the China State Administration of Work Safety issued Order No. 80, titled "Decision of the State Administration of Work Safety on Abolishing and Amending Ten Regulations in the Fields of Labor Protection Articles and Safety Training." This order explicitly repealed the previously established "Regulations on the Supervision and Administration of Labor Protection Articles" (State Administration of Work Safety Order No. 1 issued on July 22, 2005).

Article 6 of the abolished "Regulations on the Supervision and Administration of Labor Protection Products" stipulated that special labor protection products were subject to mandatory safety sign management. Consequently, the certification process for "labor safety certification" underwent a significant transformation. It shifted from the former national compulsory certification to a system of enterprise self-declaration certification.

LA-APPROVED RESPIRATORS

GB 2626-2019

Introducing the LA-approved KN95 product, featuring an ergonomic design that provides a good fit and unparalleled comfort for wearers. It offers a 95% standard minimum filtration efficiency of 0.3-micron particles, ensuring reliable breathing protection against oil-free aerosols. The adjustable nosepiece reduces the potential for eyewear fogging while providing a custom and secure seal. The inner urethane foam for the nose, combined with the high elasticity head straps, ensures a tight and secure seal for superior protection

Against non-oil-based particulates, such as minerals, iron, ore, flour or other substances.



Do not use in areas where the oxygen concentration is lower than 19.5%.
Do not use for protection against toxic gases and vapors.

LA-APPROVED RESPIRATORS

We also have dual-certified products with NIOSH N95/P95 certification



9500V-N95



9600-N95



9600V-N95



SEKURA-N95



SEKURA-N95V



MAK699+



MK710-P95VOV



MK910-N95



MKN95V-530W



2201



N95VOV-510W



6001 KN95



6001B KN95



9600V-KN95



9600V-KN95OV



6009 KN95

ITEM

DESCRIPTION

PACKAGING

6001 KN95	KN95 earloop foldable respirator(white)	1 mask/bag, 30 or 50bags/box, 12 boxes/carton
6001B KN95	KN95 earloop foldable respirator(black)	1 mask/bag, 30 or 50bags/box, 12 boxes/carton
6002 KN95	KN95 headstrap foldable respirator(white)	1 mask/bag, 30 or 50bags/box, 12 boxes/carton
6002B KN95	KN95 headstrap foldable respirator(black)	1 mask/bag, 30 or 50bags/box, 12 boxes/carton
9600-KN95	KN95 cone respirator	20 masks//box, 12 boxes/carton
9600V-KN95	KN95 cone respirator with valve	10 masks//box, 12 boxes/carton
9600-KN95OV	KN95 cone respirator for nuisance organic vapors w/valve	10 or 20 masks//box, 12 boxes/carton
6009 KN95	KN95 earloop foldable respirator(steel seal)	1 mask/bag, 30 bags/box, 12 boxes/carton



Surgical Masks

Surgical masks, widely used in various fields, including medical, dental, food processing, and laboratories, serve as examples of personal protective equipment designed to safeguard the wearer from particles or liquid contamination on the face.

Apart from offering general surgical masks, we also offer a professional medical series. These masks are designed with an additional layer of special fabric, providing effective protection against saliva and blood without compromising breathing comfort. Furthermore, these masks have successfully passed the blood penetration test of ASTM/F1862-2000a, meeting the performance level requirements of EN 14683 Type IIR standards.

NIOSH-APPROVED N95 RESPIRATORS

SURGICAL N95 RESPIRATORS

Surgical N95 respirators are examples of personal protective equipment that are used to protect the wearer from particles or from liquid contaminating the face. In addition to meeting NIOSH N95 standards, these products also pass additional testing for fluid resistance, flammability and biocompatibility. The Food and Drug Administration (FDA), Centers for Disease Control and Prevention (CDC) and the National Institute for Occupational Safety and Health (NIOSH) together regulate these respirators.




ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
9500-N95	Surgical N95 Respirator	M/L	20 masks/box; 12 boxes/case	TC-84A-5411
9500-N95S	Surgical N95 Respirator	S	20 boxes/case; 12 boxes/case	TC-84A-5463
910-N95FMX	Surgical N95 Respirator with flat fold	M/L	20 boxes/case; 12 boxes/case	TC-84A-4541
SEKURA-X	Surgical N95 Respirator with vertical fold	M/L	40 boxes/case; 12 boxes/case	TC-84A-9316

NIOSH-APPROVED RESPIRATORS

SURGICAL N95 RESPIRATORS

Brand new launches! Makrite latest surgical N95 respirators are examples of personal protective equipment that are used to protect the wearers from particles or from liquid contaminating the face. In addition to meeting NIOSH N95, CE EN149 FFP2 NR D and CE EN14683 Type IIR standards, these products also pass additional testing for fluid resistance, flammability and biocompatibility. The Food and Drug Administration (FDA) and the National Institute for Occupational Safety and Health (NIOSH) together regulate these respirators.



ITEM	DESCRIPTION	SIZE	PACKAGING	APPROVAL
9603	N95 Respirator	M/L	20 boxes/case;12 boxes/case	TC-84A-9462
9603S	N95 Respirator	S	20 boxes/case;12 boxes/case	TC-84A-9459
FIT+ 4803	N95 Respirator	M/L	20 boxes/case;12 boxes/case	TC-84A-9460

FDA 510K PLEATED FACE MASKS

Pleated mask is widely used in different industries, including medical, dental, food service, laboratories and even daily civilian uses. Here we have the same high quality pleated mask with tie-on bands and ear loops to satisfy different users' needs. Those pleated masks have passed ASTM F2100 Level 3.



CE-APPROVED & UKCA-APPROVED RESPIRATORS

CE Medical Face Masks

Medical Comfort Series are industrial and medical respirators that have passed the dual certification of European CE EN149 FFP2 NR D and CE EN14683 Type IIR, which can meet the needs of customers in both industrial and medical industry. The headstrap design is more convenient for users to adjust according to their needs. Compared with ear-loop masks, it is less likely to fall off, ensure the wearer can wear it more firmly during work, and also allow the wearer to hand it easily on their neck when not using.



9603



9603S



FIT+ 4803

ITEM	DESCRIPTION	SIZE	PACKAGING	DESCRIPTION			SIZE	PACKAGING	APPROVAL
				CE EN149	CE EN14683	Splash Resistance			
M663BE	Ear Loop Mask, Blue (ASTM F2100-19 Level 3)	M/L	50 masks/box; 40 boxes/case	FFP2 NR D	Type IIR	Pass at 21 kPa	M/L	20 boxes/case;12 boxes/case	CE 711366 & UKCA 759308
M663BT	Tie-On Mask, Blue (ASTM F2100-19 Level 3)	M/L	50 masks/box; 40 boxes/case	FFP2 NR D	Type IIR	Pass at 21 kPa	S	20 boxes/case;12 boxes/case	CE 711366 & UKCA 759308
				FFP2 NR D	Type IIR	Pass at 21 kPa	M/L	20 boxes/case;12 boxes/case	CE 711366 & UKCA 759308
				-	Type I	-	M/L	50 masks/bag; 40 bags/case	-
				-	Type II	-	M/L	50 masks/bag; 40 bags/case	-
				-	Type IIR	Pass at 16 kPa	M/L	50 masks/bag; 40 bags/case	-

ITEM	DESCRIPTION			SIZE	PACKAGING	APPROVAL
	CE EN149	CE EN14683	Splash Resistance			
9603	FFP2 NR D	Type IIR	Pass at 21 kPa	M/L	20 boxes/case;12 boxes/case	CE 711366 & UKCA 759308
9603S	FFP2 NR D	Type IIR	Pass at 21 kPa	S	20 boxes/case;12 boxes/case	CE 711366 & UKCA 759308
FIT+ 4803	FFP2 NR D	Type IIR	Pass at 21 kPa	M/L	20 boxes/case;12 boxes/case	CE 711366 & UKCA 759308
M613BE	-	Type I	-	M/L	50 masks/bag; 40 bags/case	-
M623BE	-	Type II	-	M/L	50 masks/bag; 40 bags/case	-
M633BE	-	Type IIR	Pass at 16 kPa	M/L	50 masks/bag; 40 bags/case	-



Makrite Medical Series

Disposable Medical Masks: Meet the requirements of YY/T0969-2013 Disposable Medical Masks, medical masks are suitable for wearers who are working in an ordinary medical environment with no risk of body fluids and splashes, and provides a certain physical barrier to prevent the direct penetration of pathogenic microorganisms and particulate matter.

Medical Protective Masks: Meet the requirements of GB 19083-2010 Technical Requirements for Medical Protective Masks, these masks are mainly suitable for medical institutions and people who are exposed to viral materials to prevent the transmission of viruses from patients to medical workers.

DISPOSABLE MEDICAL MASKS

Disposable medical masks are widely used in different industries, including medical treatment, food processing, electronics, pharmaceuticals, laboratories, daily life and other industries that pay high attention to hygiene. The bacterial filtration effect can reach more than 95%. This product provides health care for the wearers who work in an ordinary medical environment without risk of body fluids and splashes, creating a certain physical barrier to prevent the direct penetration of pathogenic microorganisms and particulate matter.



latex-free elastic ear loops,
allergy friendly

pleated design
increases the protection
area and reduce
pressure on the face

special static filter layer can
effectively filter dust and
microorganisms

flexible nose clip
provides excellent fit for
the majority of people

breathing resistance is
only 20 pa/cm2 for easy
breathing

skin-friendly polypropylene
non-woven fabric ensures
comfort and breathability



MEDICAL PROTECTIVE MASKS

Medical protective masks are widely used in different industries, including medical treatment, food processing, electronics, pharmaceuticals, laboratories, daily life and other industries that pay attention to hygiene. The filtering effect of particles above 0.3 microns can reach more than 95%. This product is suitable for wearers who are not surgical staff at first-line nor who are exposed to viral materials in medical institutions, which are used to prevent the spread of viruses from patients to medical workers. This series of masks have obtained both GB 19083 and NIOSH N95 certifications.



adjustable concealed
aluminum nose clip and 3D
nose foam for proper fit

breathing resistance is
only 10mmH2O for easy
breathing

special static filter layer can
effectively filter dust and
microorganisms

skin-friendly polypropylene
non-woven fabric ensures
comfort and breathability

approved by NIOSH N95

fluid-resistant at a pressure of 80
mmHg, surface moisture resistance
meets the national standard.



adjustable silicone tab on each strap
for secure personal fit



801-N95



SEKURA-X

ITEM	DESCRIPTION	SIZE	PACKAGING
M613BE	Non-sterile Ear Loop Mask (EN 14683 Type I)	17.5cm X 9.0cm	10 masks/sealed bag; 5 sealed bags/box; 12 boxes/case
M623BE	Non-sterile Ear Loop Mask (EN 14683 Type II)	17.5cm X 9.5cm	10 masks/sealed bag; 5 sealed bags/box; 12 boxes/case
M633BE	Non-sterile Ear Loop Mask (EN 14683 Type II R)	17.5cm X 9.0cm	10 masks/sealed bag; 5 sealed bags/box; 12 boxes/case
M663BE	Non-sterile Ear Loop Mask (ASTM Level III)	17.5cm X 9.5cm	10 masks/sealed bag; 5 sealed bags/box; 12 boxes/case
M663BT	Non-sterile Lace-Up Mask (ASTM Level III)	17.5cm X 9.5cm	10 masks/sealed bag; 5 sealed bags/box; 12 boxes/case

ITEM	SIZE	PACKAGING
801-N95 General Grade Vertical Folded Respirator	156mm X 120mm-4p	1 mask/sealed bag; 10 sealed bags/box; 12 boxes/case
SEKURA-X General Grade Vertical Folded Respirator	150mm x 112mm-3P	1 mask/bag; 30 bags/box; 12 boxes/case

MEDICAL PROTECTIVE MASKS

Medical protective masks are widely used in different industries, including medical treatment, food processing, electronics, pharmaceuticals, laboratories, daily life and other industries that pay attention to hygiene. The filtering effect of particles above 0.3 microns can reach more than 95%. This product is suitable for wearers who are not surgical staff at first-line nor who are exposed to viral materials in medical institutions, which are used to prevent the spread of viruses from patients to medical workers. This series of masks have obtained both GB 19083 and NIOSH N95 certifications.



ITEM	SIZE	PACKAGING
9601-N95 General Grade Arched Respirator	146mm X 117mm-5p	1 mask/sealed bag; 5 sealed bags/box; 6 boxes/ case



DISPOSABLE DUST MASKS

Mainly used for cleaning, cone-shaped dust masks are easy to put on and are shaped to cup around the nose and mouth for comfortable breathing. A malleable nosepiece adjusts for a customized fit and rounded edges help ensure a protective and comfortable fit. In addition to medical, dental and individual use, these face masks are also helpful for a variety of other applications. For example, face masks are widely used in industrial settings, the food service industry, agriculture, laboratories and DIY home improvement, among other markets. These products offer our industrial customers an economical way of expanding sales into other markets.

individually adjustable aluminum nose clip for excellent fit

ultrasonic welded headstraps for secure seal



CONE MASK FOR GENERAL USE



103B



102P



101W



DISPOSABLE PLEATED MASKS

ultrasonic welded latex-free headstraps

anti-fog transparent shield, offers eye protection from splashes

Soft PP non-woven material with low breathing resistance for comfort



Ear Loop Mask with Extended Eyeshield P801WE



Ear Loop Mask with Anti-fog Foam Strip P603BEF/P803BEF



Ear Loop Mask P603BE/P803BE



Tie-On Mask P603BT/P803BT

ITEM	DESCRIPTION	SIZE	PACKAGING
103B	Premium Dust Mask, Blue	M/L	50 masks/box; 12 or 20 boxes/case
102P	Premium Dust Mask, Pink	S	50 masks/box; 12 or 20 boxes/case
101W	Premium Dust Mask, White	M/L	50 masks/box; 12 or 20 boxes/case
203B	Premium Dust Mask, Blue	M/L	50 masks/box; 12 or 20 boxes/case
202P	Premium Dust Mask, Pink	S	50 masks/box; 12 or 20 boxes/case

ITEM	DESCRIPTION	SIZE	PACKAGING
P603BEF	Ear Loop Mask, Blue	M/L	50 masks/box; 12 or 20 boxes/case
P803BEF	Ear Loop Mask, Blue (Prevent Blood Penetration)	M/L	50 masks/box; 12 or 20 boxes/case
P602PTF	Ear Loop Mask, Pink	S	50 masks/box; 12 or 20 boxes/case
P802PTF	Ear Loop Mask, Pink (Prevent Blood Penetration)	S	50 masks/box; 12 or 20 boxes/case
P603BE	Ear Loop Mask, Blue	M/L	50 masks/box; 12 or 20 boxes/case
P803BE	Ear Loop Mask, Blue (Prevent Blood Penetration)	M/L	50 masks/box; 12 or 20 boxes/case
P602PE	Ear Loop Mask, Pink	S	50 masks/box; 12 or 20 boxes/case
P802PE	Ear Loop Mask, Pink (Prevent Blood Penetration)	S	50 masks/box; 12 or 20 boxes/case
P603BT	Tie-On Mask, Blue	M/L	50 masks/box; 12 or 20 boxes/case
P803BT	Tie-On Mask, Blue (Prevent Blood Penetration)	M/L	50 masks/box; 12 or 20 boxes/case
P602PT	Tie-On Mask, Pink	S	50 masks/box; 12 or 20 boxes/case
P802PT	Tie-On Mask, Pink (Prevent Blood Penetration)	S	50 masks/box; 12 or 20 boxes/case
P801WE	Ear Loop Mask, Blue (Prevent Blood Penetration)	M/L	50 masks/box; 12 or 20 boxes/case



Makrite & huSHE Series

Makrite is a civilian mask brand under Makrite Industries that integrates R&D, production and sales. With a mature supply chain, professional R&D strength and innovative design capabilities, it provides consumers with trendy and comfortable masks. Its product line is mainly divided into disposable face masks, stylish masks, children masks and other series.

huSHE is a civilian brand launched by Makrite in 2022. Makrite expands its expertise in manufacturing industrial masks that has been cultivated for many years to the civilian market, exploring the functional aesthetics of masks. From sports style to fashion style, consumers can always choose the one that fit their needs from a variety of products.

Test Report



报告编号: GW1381-2022

共 2 页 第 1 页

产品名称 Product name	儿童口罩	规格型号 Specification	880S
		商 标 Brand	Makrite
委托单位/地址/联系电话 Client/Addr/Tel	东莞迅安塑胶纤维制品有限公司/广东省东莞市塘厦镇 168 工业区广源街 18 号 /13431227703		
生产单位/地址/联系电话 Manufacturer/Addr/Tel	东莞迅安塑胶纤维制品有限公司/广东省东莞市塘厦镇 168 工业区广源街 18 号 /—		
样品等级 Sample grade	—	样品编号 Sample number	GW1381-2022
样品数量 Sample quantity	50 只	样品接收日期 Receiving date of sample	2022/06/22
检验检测类别 Test type	委托检验	货号/批号/款号 Article number/Batch number/Style number	062022
检验检测日期 Test date	2022/06/22 ~ 2022/06/28	检验检测地点 Test site	本公司检验室
样品状态 Sample state	符合检验检测要求	样品描述 Sample description	详见样品照片
检验检测依据 Test standard(s)	GB/T 38880-2020 《儿童口罩技术规范》		
检验检测项目 Test item(s)	甲醛含量、pH 值、颗粒物过滤效率、耐干摩擦色牢度		
检验检测结果 Test result	具体检验结果详见第 2-3 页。		
备 注 Note	样品信息由委托方提供，本报告仅对来样负责。 委托方要求所检项目按照儿童卫生口罩判定。		

批准: 苏伟群

审核: 伍慧芬

编制: 万恒



huSHE—Based on the principles of reassuring, comfort and fashion, huSHE will guard your health by improving your breathing experience.

Disposable Face Masks for Adult

- Bacterial filtration efficiency up to 95%
- Pressure-less bubble ear straps for extra comfort
- Premium melt-blown filtration layer provides upgraded protection
- Waterproof outer layer can isolate the dust and droplet better
- Skin-friendly hydrophilic inner layer ensures breathability
- Individually packaged for more convenient and easy-to-carry
- Choose from classic black and white or fashionable morandi color



Korean Style 4D Mask

- 4D design fit perfectly to the face
- Bacterial filtration efficiency up to 95%
- Pressure-less bubble ear straps for extra comfort
- Premium melt-blown filtration layer provides upgraded protection
- Waterproof outer layer can isolate the dust and droplet better
- Skin-friendly hydrophilic inner layer ensures breathability
- Individually packaged for more convenient and easy-to-carry
- Choose from classic black and white or fashionable morandi color



CF800

- Patented dual-purpose double-valves design help rapidly exhaust hot humid air, increasing breathing comfort
- Upgraded environmental concept, replaceable filter element same as N95 grade
- Washed by water without deformation
- Adjustable head strips for proper fit
- 3D shape ergonomic stylish design
- Combination of nose clip and M-shaped nose foam perfectly improves fitness of mask, reduces potential for eyewear fogging



Makrite combines product manufacturing experience, high-standard inspection and stylish aesthetic design to provide safe and high-quality masks, optimizing the wearer’s breathing experience. Classic black and white or fashionable morandi colors both are provides for easy-matching.

880S Kids Mask

This kids mask fully complies with the requirements of the China New National Standard for Children on protection performance, comfort performance and safety performance, which can meet all the needs of the children. The 3D three-dimensional structure with light-weight and breathability, leaving more space for children’s mouth and nose. The high-elastic ear loops and adjustable nose clip to ensure overall safety protection and no strain of ears. Different colorful printed design fits both boys and girls to match daily dress, which can be used indoor, school, kindergarten and other public place.

- Features:
1. The shell shape and folded design can fit different faces, wrap the chin, the gap is smaller, perfect for kids to wear.
 2. The high-elastic ear loops and adjustable nose clip to ensure overall safety protection and no strain of ears.
 3. Nosefoam-free design increase safety.
 4. Passes U.S.A Nelson authoritative sensitivity test, the skin-friendly and soft material prevent allergies.



Disposable Kids Mask

Children’s disposable masks fit comfortably and designed to be lightweight and breathable for the kids. With the feature of unfolding and expanding, the face masks will ensure full coverage of kid's nose and mouth, while the soft ear loops will not cause discomfort or irritation.

- Features:
1. Design for smaller faces provides a better seal.
 2. Three layers of filtration, comfortable and easy breathing.
 3. Extra-elastic earloops ensure a comfortable protection.
 4. Adjustable nose clip.



According to the survey of consumer demand, our R&D team developed the 9900 series which is specially designed for daily use with both professional protection and wearing comfort. Double wearing mode of earloop and headstraps allows consumers to choose the wearing method in different scenarios.



The earloop mode increase the fit and tightness by adjusting the elastic band while the headband mode liberates the ears which is not easy to fall off in the working environment. In addition to the foldable design for easy pocketing, there are also activated carbon and air valve styles for consumers to choose from.




KN95 cone mask
9900 Medium/Large



KN95 with valve
9900V Medium/Large



KN95 carbon layer for nuisance
organic vapors
9900OV Medium/Large



KN95 with valve
carbon layer for nuisance organic vapors
9900VOVMedium/Large

NEW PRODUCT LAUNCH

ITEM	DESCRIPTION	SIZE	PACKAGING
9900	KN95 Respirator	M/L	15 boxes/case;12 boxes/case
9900V	KN95 Respirator with valve	M/L	10 boxes/case;12 boxes/case
9900OV	KN95 Respirator	M/L	15 boxes/case;12 boxes/case
9900VOV	KN95 Respirator with valve	M/L	10 boxes/case;12 boxes/case

G520 NON-POWERED AIR-PURIFYING HALF MASK
PARTICULATE/CHEMICAL CARTRIDGE RESPIRATOR

Makrite G520 series of non-powered air-purifying respirators/chemical cartridge respirators have passed the inspection of the third-party institution. The respiratory mask is available against certain chemical gas or particulate when it is assembled with different Makrite’s chemical cartridges or filters, which are comply with China’s GB 2626-2019 & GB 2626-2009 non-powered air-purifying particulate respirator specifications.



The G520 series of respirators adopts a direct-connect half mask design, which is only suitable for use in the working environment that is non-air supplied. G520 series half mask respirators are available in three sizes, please refer to the product information below for the applicable head size range.

ITEM	DESCRIPTION	SIZE	APPLICABLE HEAD SIZE	PACKAGING
G520-S	Air-Purifying Half Mask Respirator	S	Face Width 120.5-144.5(MM) Face Length 98.5-118.5(MM)	1 masks/box; 24 boxes/case
G520-M	Air-Purifying Half Mask Respirator	M	Face Width 132.5-158.5(MM) Face Length 108.5-138.5(MM)	1 masks/box; 24 boxes/case
G520-L	Air-Purifying Half Mask Respirator	L	Face Width 144.5-158.5(MM) Face Length 118.5-138.5(MM)	1 masks/box; 24 boxes/case

ITEM	DESCRIPTION	PACKAGING
MF020 KN100	KN100 Particulate Filter	2pcs/bag;10bags/box;12 boxes/case
MF030 KP100	KP100 Particulate Filter	2pcs/bag;10bags/box;12 boxes/case
MWP020A KN100	Double Attachments Spray Proof Filter Retainer	1 mask/bag;10 bags/box;24 boxes/case;240pcs/case
MWP030A KP100	Double Attachments Spray Proof Filter Retainer	1 mask/bag;10 bags/box;20 boxes/case;200pcs/case
MK 9013-GB	P-A-1 Chemical Cartridge Protect Against Organic Vapors	2pcs/bag; 6bags/box; 8 boxes/case
MK 9014-GB	P-K-1 Chemical Cartridge Protect Against Ammonia and Its Derivatives	2pcs/bag; 6bags/box; 8 boxes/case
MK 9017-GB	P-E-1 Chemical Cartridge Protect Against Acid Gas	2pcs/bag;6bags/box; 8 boxes/case
MF022A KN95	Particulate Filters	2pcs/bag;50bags/box;12 boxes/case
MK9000C	Cartridge - Filter Cotton Cover	2pcs/bag;10bags/box;12 boxes/case



NEW PRODUCT LAUNCH



MK1002 COLLAPSIBLE NON-POWERED AIR-PURIFYING HALF MASK PARTICULATE RESPIRATOR

Makrite MK1002 series of non-powered air-purifying particulate respirators have passed the inspection of the third-party institution. The half mask respirator is available against particulate when it used with approved particulate filters. Comply with GB/T 18664-2002 and the U.S.A. Occupational Safety and Health Administration (OSHA) federal regulations, the respirators can be used to protect against certain airborne contaminants depending on the environment and type of respiratory hazards that need to be filtered and the respirators also comply with China’s GB 2626-2019 regulations.

ergonomically designed to perfectly fit the contours of the face


adjustable head straps are suitable for most head sizes, designed for tight and comfortable wearing

foldable silicone mask body takes up less space for easy storage

matched with particulate filter and spray proof filter retainer



This series half masks respirators has only one size. Please refer to the product information below for choosing the applicable face size. There are five color options for the body, and different colors will not affect the function and protection of the respirators.




NEW PRODUCT LAUNCH

ITEM	DESCRIPTION	COLOR	APPLICABLE HEAD SIZE	PACKAGING
MK1002A	Collapsible Air-Purifying Particulate Respirator	Black		1 masks/box; 24 boxes/case
MK1002B	Collapsible Air-Purifying Particulate Respirator	Red		1 masks/box; 24 boxes/case
MK1002C	Collapsible Air-Purifying Particulate Respirator	Blue	Face Width 132.5-158.5(MM)	1 masks/box; 24 boxes/case
MK1002D	Collapsible Air-Purifying Particulate Respirator	Green	Face Length 108.5-138.5(MM)	1 masks/box; 24 boxes/case
MK1002E	Collapsible Air-Purifying Particulate Respirator	Grey		1 masks/box; 24 boxes/case
MK1002F	Single Spray Proof Filter Retainer(matched with particulate filter)			1 pc/bag; 10 bags/box; 12 boxes/case
MF022B	KN95 Particulate Filter (can effectively protect against non-oil particle)			25 pcs/bag; 4 bags/box; 12 boxes/case; 1200 pcs/case

M200 NON-POWERED AIR-PURIFYING HALF MASK
PARTICULATE/CHEMICAL CARTRIDGE RESPIRATOR

The M200 non-powered air-purifying half mask particulate/chemical cartridge respirator series are made of advanced composite materials with high strength and high resilience. Also, the soft, lightweight material and the well-balanced design enhances comfort and allows easy cleaning. The adjustable head strap design is convenient for users in compatible with other personal protection equipment such as safety glasses, ear plugs, hard hats, etc.



The M200 half mask body can be combined with different filter accessories to meet the diverse needs of consumers.

Combination 1: combination of MK series chemical cartridges, MF022A KN95 particulate filter cotton and MK9000C filter cotton cover for protecting against certain chemical vapor and particulates.

Combination 2: combination of MK1002F single spray proof filter retainer and MF022B KN95 particulate filter cotton for protecting against particulates.

Combination 3: combination of MK1002G single filter retainer and MF022B KN95 particulate filter cotton for protecting against particulates.

ITEM	DESCRIPTION	PACKAGING
M200	Non-Powered Air-Purifying Half Mask Respirator	1 pc/bag/box; 24 boxes/case
MK2113-GB	Single P-A-1 Chemical Cartridge, Protect Against Organic Vapors (complies with GB 2890-2009 standard)	1 pc/bag; 8bags/box;10 boxes/case
MK2114-GB	Single P-K-1 Chemical Cartridge, Protect Against Ammonia and Its Derivatives (complies with GB 2890-2009 standard)	1 pc/bag; 8bags/box;20 boxes/case
MK2117-GB	Single P-E-1 Chemical Cartridge, Protect Against Acid Gas (complies with GB 2890-2009 standard)	1 pc/bag; 8bags/box;20 boxes/case
MK1002G	Single Filter Retainer (used with filter cotton)	1 pc/bag; 10bags/box;12 boxes/case
MK9000C	Cartridge - Filter Cotton Cover (used with chemical cartridge and filter cotton)	1 pc/bag; 10bags/box;12 boxes/case
MF022A	KN95 Particulate Filter Cotton (complies with GB 2890-2009 standard, used with chemical cartridge and filter cotton cover)	2 pcs/bag; 50bags/box;12 boxes/case
MF022B	KN95 Particulate Filter Cotton (complies with GB 2890-2009 standard, used with single filter cotton container)	25 pcs/bag; 4bags/box;12 boxes/case

APPLICABLE HEAD SIZE RANGE

Face Width 132.5-158.5(MM)
Face Length 108.5-138.5(MM)



NEW PRODUCT LAUNCH

G510 NON-POWERED AIR-PURIFYING HALF MASK PARTICULATE/CHEMICAL CARTRIDGE RESPIRATOR

Makrite G510 series of non-powered air-purifying particulate/chemical cartridge respirators have passed the inspection of the third-party institution. The half mask respirator is available against particulate when it used with approved particulate filters. Comply with GB/T 18664-2002 and the U.S.A. Occupational Safety and Health Administration (OSHA) federal regulations, the respirators can be used to protect against certain airborne contaminants depending on the environment and type of respiratory hazards that need to be filtered. The G510 series have approved by China’s GB 2626-2019 and NIOSH (The National Institute for Occupation Safety and Heath) regulations.

The G510 Series Non-Powered Air-Purifying Half Facepiece Respirator is made of silicon, which is less likely to cause skin swelling and itching or other allergic reactions than rubber mask.

G510 is a kind of air-purifying respirator which adopts the bayonet connection desian, only applicable for the working environment that non-powered air-purifying respirators are used.



The body of this series of half facepiece respirator adopts the design of slide-rail buckle that has obtained a US invention patent. It includes a front cover that can be assembled to or disassembled from the mask body through a linear one-way movement and the whole set of respiratory equipment is easy to disassemble and clean.

This series of half faceplece respirator adopts butterfly-shaped exhalation valve design that has obtained a US invention patent. Its performance of exhalation resistance is very well. When using, the breathing is smooth and easy to avoid the chamber becomes too hot and humid.

MK1001 COLLAPSIBLE NON-POWERED AIR-PURIFYING HALF MASK PARTICULATE RESPIRATOR

Makrite MK1001 series of non-powered air-purifying particulate respirators have passed the inspection of the third-party institution. The half mask respirator is available against particulate when it used with approved particulate filters. Comply with GB/T 18664-2002 and the U.S.A. Occupational Safety and Health Administration (OSHA) federal regulations, the respirators can be used to protect against certain airborne contaminants depending on the environment and type of respiratory hazards that need to be filtered and the respirators also comply with China’s GB 2626-2019 regulations.

free-sliding design of the head strap is convenient for the wearer to use



foldable silicone mask body takes up less space for easy storage, silicon face seal to eliminate skin swelling and itching or other allergic reactions

can match with different level particulate filter pads

The MK1001/MK1002 series of respirators adopts a direct-connect half mask design, which is only suitable for use in the working environment that is non-air supplied. The respirators are available in five colors and different colors will not affect the function and protection of the respirators, please refer to the product information below for the applicable head size range.

ITEM	DESCRIPTION	SIZE	APPLICABLE HEAD SIZE	PACKAGING
G510-S	Non-Powered Air-Purifying Particulate Respirator	S	Face Width 120.5-144.5(MM) Face Length 98.5-118.5(MM)	1 masks/box; 24 boxes/case
G510-M	Non-Powered Air-Purifying Particulate Respirator	M	Face Width 132.5-158.5(MM) Face Length 108.5-138.5(MM)	1 masks/box; 24 boxes/case
G510-L	Non-Powered Air-Purifying Particulate Respirator	L	Face Width 144.5-158.5(MM) Face Length 118.5-138.5(MM)	1 masks/box; 24 boxes/case
MF020 KN100	KN100 Particulate Filter			2pcs/bag;10bags/box;12 boxes/case
MF030 KP100	KP100 Particulate Filter			2pcs/bag;10bags/box;12 boxes/case
MK 9013-GB	P-A-1 Chemical Cartridge	Protect Against Organic Vapors		2pcs/bag; 6bags/box; 8 boxes/case
MK 9014-GB	P-K-1 Chemical Cartridge	Protect Against Ammonia and Its Derivatives		2pcs/bag; 6bags/box; 8 boxes/case
MK 9017-GB	P-E-1 Chemical Cartridge	Protect Against Acid Gas		2pcs/bag; 6bags/box; 8 boxes/case
*MF002A N95	Particulate Filters			2pcs/bag;10bags/box
MF022A KN95	Particulate Filters			2pcs/bag;10bags/box
*MK9000C	G510 - Filter Cotton Cover (used with filter cotton)			
ITEM	DESCRIPTION			PACKAGING
*MK9013	Chemical Cartridge	Protect Against Organic		2pcs/bag; 6bags/box; 8 boxes/case
*MK9017	Chemical Cartridge	Protect Against Acid Gas		2pcs/bag; 6bags/box; 8 boxes/case
*MF100	N100 Particulate Filter			2pcs/bag;10bags/box;12 boxes/case
*MF210	P100 Particulate Filter			2pcs/bag;10bags/box;12 boxes/case

*Products that have approved by NIOSH



ITEM	DESCRIPTION	COLOR	APPLICABLE HEAD SIZE	PACKAGING
MK1001A	Collapsible Air-Purifying Particulate Respirator	Black		1 masks/box; 24 boxes/case
MK1001B	Collapsible Air-Purifying Particulate Respirator	Red		1 masks/box; 24 boxes/case
MK1001C	Collapsible Air-Purifying Particulate Respirator	Blue	Face Width 132.5-158.5(MM)	1 masks/box; 24 boxes/case
MK1001D	Collapsible Air-Purifying Particulate Respirator	Green	Face Length 108.5-138.5(MM)	1 masks/box; 24 boxes/case
MK1001E	Collapsible Air-Purifying Particulate Respirator	Grey		1 masks/box; 24 boxes/case
MF020	KN100 Particulate Filter			2pcs/bag; 10bags/box; 12 boxes/case
MF030	KP100 Particulate Filter			2pcs/bag; 10bags/box; 12 boxes/case

WHY CHOOSE MAKRITE?

As with any business decision, there are plenty of factors to consider. We see these as among the most important:



Highest-quality products

Our top priority is to provide you with the best possible disposable respirators. We do this by investing a lot of time and energy in our initial product designs and development. We then ensure that each and every respirator is produced to our exacting standards using sophisticated, specially designed equipment. Our quality control and laboratory testing departments complete the process so that we consistently ship the same product every time you order.



Globally competitive pricing

Our pricing is value-oriented, guaranteeing a high-quality product that meets rigorous NIOSH and other government-approved standards at very competitive prices in the US and elsewhere around the world. And while we may not be the lowest priced manufacturer in Asia, given the number of major companies that turn to us for private brand respirators, we are the most significant and reliable one.



Private brand packaging

Many of our customers market our products in their own brand packaging. We encourage this since it enables you to sell the best respirators in the world with your company's name on the mask and packaging. When you order from Makrite, both the smaller inner boxes with the masks and the outer shipping cartons are branded with your company's name. Plus the only place your customers can purchase the product is through your company.



Unbeatable customer service

At Makrite we offer only the highest levels of service and support to customers around the globe. We provide prompt and personal responses to all of your questions, comments and concerns. Regardless of where you are or how large or small your concern might be, we invite you to call us at **1-800-379-9929** or email us at **info@makritenorthamerica.com**. We take care of the rest.



Unmatched technical support

When you do business with Makrite, we are your support team. We work diligently to answer all of your technical questions within 24 hours or less.




Excellent marketing & sales support

We also help your business by offering services such as product training programs, product photos, catalog preparation assistance, strategic discussions and more.



Retail packaging options

All of our respirators and other face mask products are available in a wide range of packaging options to meet your needs—and those of your customers. These include clam shells, polybags with header cards, blister cards and more. We also offer customizable vending packs of one or two respirators in a bag with an instruction sheet.







Makrite Industries Inc.

TAIWAN HEAD OFFICE

11F-5, No.79, Sec. 1, Xintai 5th Rd., XiZhi Dist.,
New Taipei City 221, Taiwan
Tel: +886 2 2698 2419
Fax: +886 2 2698 2423

SHANGHAI OFFICE

Room 512, 5th Floor, Building 12, No.800, Pu Xing Rd.,
Pu Jiang Town, Min Hang District, Shanghai, China
Tel: +86 021 5431 3117
Fax: +86 021 5431 3119

HUBEI OFFICE

NO.1 Zonger Rd., Guanggu Industrial Park,
Gaoxin District, Chibi City, Hubei Province, China
Tel: +86 0715 590 0905
Fax: +86 0715 590 0906

THAILAND OFFICE

419/51 Moo. 10, Khao Khan song, Si Racha,
Chonburi 20110, Thailand
Tel: +66 033 640 868 Dial 101
Fax: +66 033 640 966



www.makrite.com



www.makrite.com.cn
sales@makrite.com